



# SLIMWAVESYSTEM

## *Slim Wave to GO System Operations Manual*



PHI Management Inc. 12121 Little Road #323 Hudson, FL 34667  
Call or Text: (813) 507-4429 Email: [Sales@SlimWaveSystem.com](mailto:Sales@SlimWaveSystem.com)



# SLIM WAVESYSTEM

**Reinvent your BODY  
Reinvent your LIFE.**



## User Instructions:

- Simply Stretch out the ends of the belt and place it comfortably on bare skin around the body area desired.
- Power on the belt using the power button on the front of the unit
- Automatic vibration enhances the effects of the LED light.
- Charge using the USB Cord provided using the USB port on the inside on the belt (*Charging requires 3 hours. Battery Life 2.5 hours of use.*)
- *Register your unit at:*  
[www.slimwavesystem.com/register-purchase/](http://www.slimwavesystem.com/register-purchase/)

**The Slim Wave to GO System uses LED lights and vibration to diminish body fat. Simply place the adjustable LED light belt around your waist on your bare skin and turn the unit on. The red LED light penetrates skin and opens the pores of fat cells allowing fatty acids to escape and burn off more easily during exercise. Use 15 minutes per treatment area, 2 to 3 times per week in addition to exercise and a healthy diet.**

The same LED cellulite reduction system used in medical offices, is now available for the home and business use. The Slim Wave to Go System uses LED lights to open up fat cells, so fatty acids can be burned off more easily during exercise.



**To get the best results follow these simple guidelines, this will enhance the effectiveness of the Slim Wave to GO System.**

- Make sure to track your food consumption continually throughout the day
- Drink Plenty of Water after using the device (Ideally you want to drink ½ ounce of water for every pound of body weight every day. For example a 130lb person should drink at least 65 ounces of WATER each day)
- Make sure to exercise for at least 30 minutes each day (approximately 500 calories of exercise)
- Do not eat 2 hours before or one hour after each session for best results

**Maintenance:** Cleaning and disinfection: After each use, using a soft cloth dampened with warm water or isopropyl alcohol, wipe the entire face of the FLH including the plastic window of the light-emitting surface. **CAUTION:** Do not allow liquid to enter any of the openings on the unit.

### **Preventative Inspection:**

- Before each use, carefully inspect the Slim Wave to GO System to make sure it does not appear damaged.
- Before operation the user should check to make sure the protective cover housing of the LED is intact.
- The Slim Wave to GO System has no serviceable parts or calibration requirements.

### **Specifications:**

#### Power Requirements

Power: 7.4v Battery / Charger: Mini USB 3.7v

Output: 10v 1a

Led: Class 1 Non-Prescription

Operating Temperature: 50° / 104 °F (10° / 40° C)

Storage Temperature: -20° / 130° F (30° -55° C)

Relative Humidity: 10- 95% storage, non condensing

Storage: Whenever possible the Slim Wave to GO System should be stored at room temperature in a cool dry place. Bathroom storage is not recommended.



### **Contraindications:**

Below is a list of contraindications. Contraindications are conditions or factors that serve as a reason not to use this machine, as it could be harmful to the individual using the SlimWaveSystem. If at any time there is a question as to whether the SlimWaveSystem should be used, consult a medical doctor before use.

Not for use on:

- Open wounds or sores
- Minors (those below 18 years of age)
- People with a pacemaker
- People who have photosensitivity
- People with Epilepsy
- People with Cancer
- People with severe kidney or liver problems
- People with thyroid problems
- Women who are pregnant or nursing

**Warranty & Service Information** Register your device at: [www.slimwavesystem.com/register-purchase/](http://www.slimwavesystem.com/register-purchase/)

### **Proprietary Notice:**

Information contained in this document is copyrighted by PHI MGMT, and may not be duplicated in full or part by any person without prior written approval of PHI MGMT. Its purpose is to provide the user with adequately detailed documentation to efficiently install, operate and maintain the Slim Wave System. Efforts have been made to keep the information contained in this document current and accurate as of the date of publication or revision.

### **Limited Warranty:**

PHI MGMT warrants the Slim Wave System to be free from defects in workmanship and material under normal use for a period of one (1) year from the date of shipment. The sole obligation of PHI MGMT under this warranty will be to repair or replace the system, at its option, should the Slim Wave System prove to be defective during the warranty period. The foregoing shall be the sole warranty remedy. As set forth herein, PHI MGMT makes no warranties, either expressed or implied, including the implied warranties of merchantability and fitness for a specific purpose. No warranty is provided if the Slim Wave System is modified without the express written consent of PHI Management. PHI Management shall not be liable in any event for incidental consequential damage. This warranty is not activated before the customer returns the completed Slim Wave System registration form within 30 days after the purchase of the product. This warranty is not assignable.

### **Service Support:**

To return your Slim Wave System for repair/replacement under the terms of this one year limited warranty, you must contact us for a return merchandise authorization number (RMA number) within one (1) year from the shipping date. Email us at [sw@slimwavesystem.com](mailto:sw@slimwavesystem.com) to obtain the RMA number from our customer service department. The return of the Slim Wave System and all accompanying documents must be packed carefully to prevent shipping damage. The RMA number must be clearly printed on the outside of the shipping container.

### **Warnings, Cautions & Notes:**

- **Warning:** The operator must be thoroughly familiar with the information in this manual before using the device.
- **Warning:** Not intended for ocular applications or direct eye exposure
- **Warning:** Use carefully. May cause serious burns.
- **Warning:** To avoid electrical shock hazard, do not open the base.
- **Warning:** Opening the base voids all warranties.
- **Warning:** Not for use on open wounds or sores, minors, people with pacemakers, photosensitivity, epilepsy, and form of cancer, kidney problems, liver problems, thyroid problems, or women who are pregnant or nursing. See CONTRAINDICATIONS section for more detailed information.
- **Warning:** The unattended use of Slim Wave System by children or incapacitated persons may be dangerous.
  
- **Caution:** Do not autotave or immerse light source in liquid.
  
- **Note:** The Slim Wave System is a Class 1 LED Device, which is safe under reasonably foreseeable conditions (per IEC 60825-1:1993+A1:1997+A2:2001).
- **Note:** The equipment is suitable for use in medical offices or households.
- **Note:** Use of mobile phones or other electronic devices around the Slim Wave Device may interfere with the device and could cause a malfunction.
- **Note:** Use or storage of this device in an environment of high humidity or precipitation may adversely affect the device.  
We advise storing in a cool dry location, not in a bathroom.
- **Note:** There are no customer serviceable parts in this device. If service is required, please contact PHI Management at [slimwavesystem@gmail.com](mailto:slimwavesystem@gmail.com) for instructions.

**Please Note: This is a personal care item and is returnable only if in its original unopened packaging and unused.**